

Anti-Lpar1 antibody

Catalog Number: 175011

Product name

Polyclonal

Anti-Lpar1 antibody

Ig Type**Specificity**

Rabbit IgG

Mouse

Applications**Antibody description**

WB, IHC-P

Rabbit polyclonal antibody to Lpar1


Dilutions**Preparation**

WB: 1:1,000

This antigen of this antibody was synthetic peptide within n-terminal mouse lpa receptor 1.

Validations**Formulation**

Liquid, 1*PBS (pH7.4), 0.2% BSA, 25% Glycerol.
Preservative: 0.05% Sodium Azide.

**Storage**

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Fig1: Western blot analysis on mouse liver using anti-LPA-receptor I polyclonal antibody.

Clonality