

Anti-FFAR2 antibody



Catalog Number: 176011

Product name

Polyclonal

Anti-FFAR2 antibody

Ig Type**Specificity**

Rabbit IgG

Human

Applications**Antibody description**

IF

Rabbit polyclonal antibody to FFAR2

Dilutions**Preparation**

IF: 1:100-1:500

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human ffar2. aa range:271-320

Validations**Formulation**

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.


Clonality

Fig1.; Immunofluorescence analysis of COS7 cells, using FFAR2 Antibody. The picture on the right is blocked with the synthesized peptide.