

Anti-FFAR2 antibody

Catalog Number: 176011

Product name

Anti-FFAR2 antibody

Specificity

Human

Antibody description

Rabbit polyclonal antibody to FFAR2

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human ffar2. aa range:271-320

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

IF

Dilutions

IF: 1:100-1:500


Validations

Fig1.; Immunofluorescence analysis of COS7 cells, using FFAR2 Antibody. The picture on the right is blocked with the synthesized peptide.