

# TFF2 antibody

Catalog Number: 116011

## Product name

TFF2 antibody

IHC: 1:20-1:200


## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

TFF2 Rabbit Polyclonal antibody. Positive IHC detected in human stomach cancer tissue, human stomach tissue, mouse liver tissue, mouse pancreas tissue. Positive WB detected in human stomach tissue. Observed molecular weight by Western-blot: 18-20 kDa

## Validations




## Preparation

This antibody was obtained by immunization of TFF2 recombinant protein (Accession Number: NM\_005423). Purification method: Antigen affinity purified.

human stomach tissue were subjected to SDS PAGE followed by western blot with Catalog No:116011(TFF2 antibody) at dilution of 1:100

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.



## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal


Immunohistochemical of paraffin-embedded human stomach cancer using Catalog No:116011(TFF2 antibody) at dilution of 1:50 (under 40x lens)

## Ig Type

Rabbit IgG

## Applications

ELISA, IHC, WB



## Dilutions

Recommended Dilution:

WB: 1:200-1:1000

Immunohistochemical of paraffin-embedded human stomach cancer using Catalog

# TFF2 antibody



Catalog Number: 116011

---

No:116011(TFF2 antibody) at dilution of 1:50 (under 10x lens)