

ZNF606 antibody

Catalog Number: 117011

Product name

ZNF606 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

ZNF606 Rabbit Polyclonal antibody. Positive IP detected in mouse liver tissue. Positive WB detected in mouse liver tissue. Observed molecular weight by Western-blot: 70 kDa

Preparation

This antibody was obtained by immunization of ZNF606 recombinant protein (Accession Number: NM_001348022). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IP, WB


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:1000

Validations



mouse liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:117011(ZNF606 Antibody) at dilution of 1:600



IP Result of anti-ZNF606 (IP:Catalog No:117011, 4ug; Detection:Catalog No:117011 1:300) with mouse liver tissue lysate 4000ug.