

CD23 antibody

Catalog Number: 109011

Product name

CD23 antibody

Specificity

Human; other species not tested.

Antibody description

CD23 Rabbit Polyclonal antibody. Positive IHC detected in human tonsillitis tissue, human lymphoma tissue.

Preparation

This antibody was obtained by immunization of CD23 recombinant protein (Accession Number: BC064417). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IHC

Dilutions


Recommended Dilution:

IHC: 1:50-1:500

Validations



Immunohistochemistry of paraffin-embedded human tonsillitis tissue slide using Catalog No:109011 (CD23, FCER2 Antibody) at dilution of 1:200 (under 10x lens). Heat-mediated antigen retrieval with Tris-EDTA buffer (pH 9).



Immunohistochemistry of paraffin-embedded human tonsillitis tissue slide using Catalog No:109011 (CD23, FCER2 Antibody) at dilution of 1:200 (under 40x lens). Heat-mediated antigen retrieval with Tris-EDTA buffer (pH 9).