

EXOSC4 antibody

Catalog Number: 110411

Product name

EXOSC4 antibody

IHC: 1:20-1:200

IF: 1:10-1:100


Specificity

Human, Mouse, Rat; other species not tested.

Validations

Antibody description

EXOSC4 Rabbit Polyclonal antibody. Positive WB detected in Jurkat cells, human heart tissue, human liver tissue, mouse liver tissue, rat liver tissue. Positive IF detected in HepG2 cells. Positive IHC detected in human skin cancer tissue. Observed molecular weight by Western-blot: 27kd




Preparation

This antibody was obtained by immunization of EXOSC4 recombinant protein (Accession Number: NM_019037). Purification method: Antigen affinity purified.

Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:110411(EXOSC4 antibody) at dilution of 1:500

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.



Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Immunohistochemical of paraffin-embedded human skin cancer using Catalog No:110411(EXOSC4 antibody) at dilution of 1:100 (under 10x lens)

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC, IF

Dilutions

Recommended Dilution:

WB: 1:500-1:5000




Immunohistochemical of paraffin-embedded

EXOSC4 antibody

Catalog Number: 110411

human skin cancer using Catalog
No:110411(EXOSC4 antibody) at dilution of 1:100
(under 40x lens)



Immunofluorescent analysis of HepG2 cells, using
EXOSC4 antibody Catalog No:110411 at 1:25
dilution and Rhodamine-labeled goat anti-rabbit
IgG (red).