

HJURP antibody



Catalog Number: 111411

Product name

HJURP antibody

Specificity

human; other species not tested.

Antibody description

HJURP Rabbit Polyclonal antibody. Positive WB detected in A549 cells, A431 cells, HeLa cells. Observed molecular weight by Western-blot: 70kd,83-85kd

Preparation

This antibody was obtained by immunization of HJURP recombinant protein (Accession Number: NM_018410). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



A549 cells were subjected to SDS PAGE followed by western blot with Catalog No:111411(HJURP antibody) at dilution of 1:800