

Anti-FGF6 antibody

Catalog Number: 176411

Product name

Anti-FGF6 antibody

Specificity

Human, Mouse, Rat

Antibody description

Rabbit polyclonal antibody to FGF6

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from the c-terminal region of human fgf6. aa range:159-208

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB

Dilutions

WB: 1:500-1:2000

Validations




Fig1.; Western Blot analysis of 293, HepG2, PC-3, M21 cells using FGF-6 Polyclonal Antibody.. Secondary antibody (catalog#: HA1001) was diluted at 1:20000