

# Anti-FGF6 antibody



Catalog Number: 176411

## Product name

Polyclonal

Anti-FGF6 antibody

## Ig Type

## Specificity

Rabbit IgG

Human, Mouse, Rat

## Applications

## Antibody description

WB

Rabbit polyclonal antibody to FGF6

## Dilutions

## Preparation

WB: 1:500-1:2000

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from the c-terminal region of human fgf6. aa range:159-208

## Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

## Clonality

## Validations




Fig1.; Western Blot analysis of 293, HepG2, PC-3, M21 cells using FGF-6 Polyclonal Antibody.. Secondary antibody (catalog#: HA1001) was diluted at 1:20000