

Anti-JAK1 (8B8) Mouse antibody



Catalog Number: 168411

Product name

Anti-JAK1 (8B8) Mouse antibody

Specificity

Human, Rat

Antibody description

JAK1 (8B8) Mouse monoclonal antibody

Preparation

Antigen: Purified recombinant human Jak1 protein fragments expressed in E.coli.

Formulation

PBS(pH 7.4) containing with 0.03% Proclin300 and 50% glycerol.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Mouse IgG2a

Applications


WB, ICC/IF, IP

Dilutions


WB: 1/1000

;ICC: 1/200

Validations




Immunocytochemistry staining of HeLa cells fixed with 4% Paraformaldehyde and using anti-Jak1



mouse mAb (dilution 1:200).

Western blot analysis of extracts from C6,Ramos,Jurkat and HeLa cell lysates using Jak1 mouse mAb (1:1000 diluted).Predicted band size:130KDa.Observed band size:130KDa.



Immunoprecipitation analysis of HeLa cell lysates using Jak1 mouse mAb.