

# INTS9 antibody

Catalog Number: 111811

## Product name

INTS9 antibody

## Specificity

Human; other species not tested.

## Antibody description

INTS9 Rabbit Polyclonal antibody. Positive WB detected in Jurkat cells, HEK-293 cells, HeLa cells. Positive IP detected in HeLa cells. Observed molecular weight by Western-blot: 74 kDa

## Preparation

This antibody was obtained by immunization of INTS9 recombinant protein (Accession Number: NM\_018250). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:500-1:5000

## Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:111811(INTS9 antibody) at dilution of 1:1000



IP Result of anti-INTS9 (IP:Catalog No:111811, 4ug; Detection:Catalog No:111811 1:1000) with HeLa cells lysate 880ug.