

# TAC3 Polyclonal Antibody



Catalog Number: 164811

## Product name

TAC3 Polyclonal Antibody

## Specificity

Human, Mouse, Rat

## Antibody description

Polyclonal antibody to TAC3

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-121 of human TAC3 (NP\_037383.1).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


WB, IF

## Dilutions

WB 1:500 - 1:2000


IF 1:10 - 1:100

## Validations



Western blot - TAC3 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using TAC3 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure



time: 20s.

Immunofluorescence - TAC3 Polyclonal Antibody

Immunofluorescence analysis of HeLa cells using TAC3 antibody. Blue: DAPI for nuclear staining.