

Anti-NMBR antibody



Catalog Number: 175811

Product name

Anti-NMBR antibody

Specificity

Human, Mouse, Rat

Antibody description

Rabbit polyclonal antibody to NMBR

Preparation

This antigen of this antibody was synthetic peptide of nmbr

Formulation

Liquid, PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P

Dilutions

WB: 1:1000-1:2000

IHC-P: 1:200-1:500

Validations




Fig1:: Western blot analysis of 1) Human Brain Tissue, 2) Rat Brain Tissue using NMBR Polyclonal Antibody.. Secondary antibody [catalog#: HA1001) was diluted at 1:20000




Fig2:: Immunohistochemical analysis of paraffin-embedded Rat Brain Tissue using NMBR Polyclonal Antibody.