

STXBP3 antibody

Catalog Number: 115811

Product name

STXBP3 antibody

IP: 1:200-1:2000

IHC: 1:20-1:200

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

STXBP3 Rabbit Polyclonal antibody. Positive IHC detected in human lymphoma tissue. Positive IP detected in K-562 cells. Positive WB detected in HepG2 cells, K-562 cells. Observed molecular weight by Western-blot: 68kd

Preparation

This antibody was obtained by immunization of STXBP3 recombinant protein (Accession Number: NM_007269). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, IHC, IP

Dilutions


Recommended Dilution:

WB: 1:500-1:5000


Validations



HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:115811(STXBP3 antibody) at dilution of 1:800



Immunohistochemical of paraffin-embedded human lymphoma using Catalog No:115811(STXBP3 antibody) at dilution of 1:100 (under 10x lens)



Immunohistochemical of paraffin-embedded

STXBP3 antibody

Catalog Number: 115811

human lymphoma using Catalog No:115811(STXBP3 antibody) at dilution of 1:100 (under 40x lens)

IP Result of anti-STXBP3 (IP:Catalog No:115811, 4ug; Detection:Catalog No:115811 1:500) with K-562 cells lysate 1200ug.

