

# Anti-PTEN Rabbit antibody



Catalog Number: 166811

## Product name

Anti-PTEN Rabbit antibody

## Specificity

Human, Mouse, Rat

## Antibody description

PTEN Rabbit polyclonal antibody

## Preparation

Antigen: A synthesized peptide derived from human PTEN

## Formulation

Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

WB, IHC-P, ICC/IF


## Dilutions

WB: 1/1000-1/2000


IHC: 1/50-1/200

ICC/IF: 1/50-1/200


## Validations



Immunohistochemical analysis of paraffin-embedded human colon, using PTEN Antibody.



Immunofluorescent analysis of HeLa cells, using



PTEN Antibody.

Western blot analysis of extracts 293T and HeLa cells using PTEN Rabbit pAb at 1:1000 dilution. Predicted band size: 54kDa. Observed band size: 54kDa.