

Recombinant Human KRAS (C118S) Protein



Catalog Number: 552111

Protein Construction

DNA encoding 1-169 aa (P01116-2, C118S) was fused with a 6-His tag.

Organism

Human

Expression Host

E.coli

Activity

Not tested.

Purity

>85%, by SDS-PAGE with Coomassie Brilliant Blue staining.

Endotoxin

Not tested.

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Met

Molecular Mass

Predictes a molecular mass of 21.53 kDa. In SDS-PAGE under reducing conditions, it migrates as an approximately 23.7 kDa band.

Formulation

Lyophilized from sterile 25 mM HEPES, 150 mM NaCl, pH 7.4


Storage

Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Reconstitute at 0.25 µg/µl in sterile water for short-term storage. Reconstitution with 50% glycerol solution is recommended for longer term storage (see Stability and Storage for more details). If a different concentration is needed for your purposes please adjust the reconstitution volume as required (please note: the ion concentration of the final solution will vary according to the volume used). Note: Centrifuge vial before opening. When reconstituting, gently pipet and wash down the sides of the vial to ensure full recovery of the protein into solution.

SDS-PAGE



Recombinant Human KRAS (C118S) Protein SDS-PAGE

M: Marker; Lane 1: KRAS (C118S); Lane 2: BSA