

Human PSMD10 recombinant protein



Catalog Number: 518111

General Information

Gene Name Synonym

26S Proteasome Non-ATPase Regulatory Subunit 10, 26S Proteasome Regulatory Subunit p28, Gankyrin, p28(GANK), PSMD10

Protein Construction

A DNA sequence encoding the Human PSMD10 Met1-Gly226 is expressed.

Organism

Human

Expression Host

E. coli

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μg protein as determined by the LAL

method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a $0.2\ \mu\text{m}$ filtered solution of 20mM PB, 150mM NaCl, 10% Glycerol, pH 7.4.

Usage Guide

Storage

Store at $<-20^{\circ}\text{C}$, minimize freeze-thaw cycles.

SDS-PAGE



Human PSMD10 recombinant protein