

RAD18 antibody

Catalog Number: 114511

Product name

RAD18 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

RAD18 Rabbit Polyclonal antibody. Positive WB detected in K-562 cells, A549 cells, HeLa cells, MCF7 cells. Observed molecular weight by Western-blot: 75kd, 85kd

Preparation

This antibody was obtained by immunization of RAD18 recombinant protein (Accession Number: NM_020165). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, IP

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



K-562 cells were subjected to SDS PAGE followed by western blot with Catalog No:114511(RAD18 antibody) at dilution of 1:300