

Phospho-EphA3-Y779 pAb



Catalog Number: 166511

Product name

Phospho-EphA3-Y779 pAb

Specificity

Human

Antibody description

Polyclonal antibody to Phospho-EphA3-Y779

Preparation

Antigen: A synthetic phosphorylated peptide around Y779 of human EphA3 (NP_005224.2).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

WB

Dilutions

WB 1:500 - 1:2000

Validations



Western blot - Phospho-EphA3-Y779 pAb

Western blot analysis of extracts of 293T and HeLa cells, using Phospho-EphA3-Y779 antibody at 1:2000 dilution. 293T cells were treated by Pervanadate (1mM) for 30 minutes after serum-starvation overnight. HeLa cells were treated by Pervanadate (1nM) for 30 minutes after serum-starvation overnight. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% BSA. Detection: ECL Basic Kit. Exposure time: 1s.