

Anti-SERPINB3 antibody



Catalog Number: 142911

Product name

Anti-SERPINB3 antibody

Immunogen

[Human SERPINB3 / SCCA-1 \(His Tag\) recombinant protein](#)

Specificity

Human SERPINB3 / SCCA-1

Antibody description

Rabbit Monoclonal to Human SERPINB3

Preparation

This antibody was obtained from a rabbit immunized with purified, recombinant Human SERPINB3 / SCCA-1 .

Formulation

0.2 µm filtered solution in PBS

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid

contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

Clonality

Monoclonal

Ig Type

Rabbit IgG


Applications

IP

Dilutions

IP: 4-8 µl/mg of lysate

Validations



Item	Lanes	A
Sample (whole cell lysate)		A431
Sample quantity		0.5 mg
IP antibody quantity		4 µl.
Protein G agarose		15 µl of 50% Protein G agarose
Gel		13% SDS-PAGE reducing gel
Primary antibody		SERPINB3 antibody at 10 µg/ml
Secondary antibody		Dylight 800-labeled antibody to rabbit IgG (F1+L), at 1:5000 dilution

Human SERPINB3 / SCCA-1

Immunoprecipitation(IP) 15631