

# AGPAT1 Polyclonal Antibody



Catalog Number: 164911

## Product name

AGPAT1 Polyclonal Antibody

## Specificity

Human, Mouse, Rat

## Antibody description

Polyclonal antibody to AGPAT1

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 204-283 of human AGPAT1 (NP\_006402.1).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


WB, IP

## Dilutions

WB 1:500 - 1:2000


IP 1:50 - 1:100

## Validations



Western blot - AGPAT1 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using AGPAT1 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL



Basic Kit .Exposure time: 90s.

Immunoprecipitation - AGPAT1 Polyclonal Antibody

Immunoprecipitation analysis of 200ug extracts of 293T cells, using 3 ug AGPAT1 antibody . Western blot was performed from the immunoprecipitate using AGPAT1 antibody at a dilution of 1:1000.