

TCL1A antibody

Catalog Number: 115911

Product name

TCL1A antibody

Specificity

Human; other species not tested.

Antibody description

TCL1A Rabbit Polyclonal antibody. Positive WB detected in Raji cells. Positive IHC detected in human tonsillitis tissue, human lung cancer tissue. Observed molecular weight by Western-blot: 14-16kd

Preparation

This antibody was obtained by immunization of TCL1A recombinant protein (Accession Number: XM_017021676). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IHC: 1:20-1:200


Validations



Raji cells were subjected to SDS PAGE followed by western blot with Catalog No:115911(TCL1A antibody) at dilution of 1:300



Immunohistochemistry of paraffin-embedded human tonsillitis tissue slide using Catalog No:115911(TCL1A Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human tonsillitis tissue slide using Catalog No:115911(TCL1A Antibody) at dilution of 1:50

TCL1A antibody



Catalog Number: 115911

(under 40x lens)