

Anti-GDF15 Rabbit antibody



Catalog Number: 166911

Product name

Anti-GDF15 Rabbit antibody

Specificity

Human

Antibody description

GDF15 Rabbit polyclonal antibody

Preparation

Antigen: A synthetic peptide of human GDF15

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

WB

Dilutions

WB: 1/1000

Validations



Western blot detection of GDF15 in HeLa cell lysates using GDF15 Rabbit pAb(1:1000 diluted).Predicted band size:34kDa.Observed band size:34kDa.