

Anti-Amyloid Precursor Protein Rabbit antibody



Catalog Number: 168911

Product name

Anti-Amyloid Precursor Protein Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

Amyloid Precursor Protein Rabbit monoclonal antibody

Preparation

Antigen: A synthetic peptide of human Amyloid Precursor Protein

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, ICC/IF, FC, IP

Dilutions

1:1000


1:20-1:500

1:20

1:20

1:20

Validations



Western blot detection of Amyloid Precursor Protein in Rat Brain, 3T3, HeLa cell lysates using Amyloid Precursor Protein Rabbit mAb (1:1000 diluted). Predicted band size: 87 kDa. Observed band size: 100 kDa.