

Anti-FCGR2 antibody



Catalog Number: 141211

Product name

Anti-FCGR2 antibody

Immunogen

[Human CD32a / Fc gamma RIIA / FCGR2A \(167 Arg, His Tag\) recombinant protein](#)

Specificity

Human CD32a

Has cross-reactivity with Human CD32b

Antibody description

Mouse Monoclonal to Human FCGR2

Preparation

This antibody was produced from a hybridoma resulting from the fusion of a mouse myeloma with B cells obtained from a mouse immunized with purified, recombinant Human CD32a/FCGR2A/Fc gamma RIIA and conjugated with APC under optimum conditions, the unreacted APC was removed.

Formulation

Aqueous solution containing 0.5% BSA and 0.09% sodium azide

Storage

This antibody is stable for 12 months from date of receipt when stored at 2°C-8°C. Protected from prolonged exposure to light. Do not freeze ! Sodium azide is toxic to cells and should be disposed of properly. Flush with large volumes of water during disposal.

Clonality

Monoclonal

Ig Type

Mouse IgG1


Applications

FCM

Dilutions

FCM

Validations



Human CD32a/FCGR2A/Fc gamma RIIA Flow Cytometry (FC) 15894