

# GRK1 antibody

Catalog Number: 111211

## Product name

GRK1 antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

GRK1 Rabbit Polyclonal antibody. Positive IP detected in Y79 cells. Positive WB detected in Y79 cells, mouse eye tissue. Observed molecular weight by Western-blot: 64 kDa

## Preparation

This antibody was obtained by immunization of GRK1 recombinant protein (Accession Number: NM\_002929). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:2000

## Validations



Y79 cells were subjected to SDS PAGE followed by western blot with Catalog No:111211(GRK1 Antibody) at dilution of 1:600



IP Result of anti-GRK1 (IP:Catalog No:111211, 4ug; Detection:Catalog No:111211 1:500) with Y79 cells lysate 1240ug.