

LGR5 antibody

Catalog Number: 112211

Product name

LGR5 antibody

Specificity

Human; other species not tested.

Antibody description

LGR5 Rabbit Polyclonal antibody. Positive IHC detected in human small intestine tissue, human colon cancer tissue, human colon tissue, human liver cancer tissue.

Preparation

This antibody was obtained by immunization of LGR5 recombinant protein (Accession Number: NM_003667). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IHC

Dilutions


Recommended Dilution:

IHC: 1:20-1:200

Validations



Immunohistochemical of paraffin-embedded human small intestine using Catalog No:112211(LGR5 antibody) at dilution of 1:200 (under 10x lens)



Immunohistochemical of paraffin-embedded human small intestine using Catalog No:112211(LGR5 antibody) at dilution of 1:200 (under 40x lens)