

# DACH1 Polyclonal Antibody



Catalog Number: 164211

---

**Product name**

DACH1 Polyclonal Antibody

**Specificity**

Human, Mouse, Rat

**Antibody description**

Polyclonal antibody to DACH1

**Preparation**

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 450-570 of human DACH1 (NP\_542937.2).

**Formulation**

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Storage**

Store at -20°C. Avoid freeze / thaw cycles.

**Clonality**

Polyclonal

**Ig Type**

Rabbit IgG


**Applications**

WB, IHC

**Dilutions**

WB 1:200 - 1:2000

IHC 1:20 - 1:200

**Validations**

Immunohistochemistry - DACH1 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human colon carcinoma using DACH1 antibody at dilution of 1:100 (40x lens).