

Mouse EPHA3/Eph Receptor A3 (His Tag) recombinant protein



Catalog Number: 504211

General Information

Protein Construction

A DNA sequence encoding the mouse EPHA3 (NP_034270.1) (Met1-His541) was expressed with a C-terminal polyhistidine tag.

Organism

Mouse

Expression Host

Human Cells

QC Testing

Activity

1. Measured by its binding ability in a functional ELISA.
2. Immobilized mouse EPHA3-His at 10 µg/mL (100 µL/well) can bind mouse EFNA5-Fc (Cat:502960). The EC₅₀ of mouse EFNA5-Fc (Cat:502960) is 4.9-11.4ng/mL.

Purity

> 95 % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per µg of the protein as determined by the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Glu 21

Molecular Mass

The recombinant mouse EPHA3 comprises 532 amino acids and has a predicted molecular mass of 60.2 kDa.

Formulation

Lyophilized from sterile PBS, pH 7.4.

1. 5 % trehalose and mannitol are added as protectants before lyophilization.
2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Mouse EphA3 Protein (His Tag) SDS-PAGE