

# PROC antibody



Catalog Number: 114211

---

## Product name

PROC antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

PROC Rabbit Polyclonal antibody. Positive WB detected in mouse liver tissue. Observed molecular weight by Western-blot: 52 kDa

## Preparation

This antibody was obtained by immunization of PR<sup>o</sup>C recombinant protein (Accession Number: NM\_001375613). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



mouse liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:114211(PRoc Antibody) at dilution of 1:1000