

ASB14 antibody



Catalog Number: 108211

Product name

ASB14 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

ASB14 Rabbit Polyclonal antibody. Positive IHC detected in human tonsillitis tissue. Positive WB detected in mouse spleen tissue, Raji cells. Observed molecular weight by Western-blot: 65 kDa

Preparation

This antibody was obtained by immunization of ASB14 recombinant protein (Accession Number: BC114953). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

Validations



mouse spleen tissue were subjected to SDS PAGE followed by western blot with Catalog No:108211(ASB14 Antibody) at dilution of 1:600



Immunohistochemistry of paraffin-embedded human tonsillitis tissue slide using Catalog No:108211(ASB14 Antibody) at dilution of 1:50 (under 40x lens)