

MICALL2 antibody



Catalog Number: 112611

Product name

MICALL2 antibody

Specificity

Human; other species not tested.

Antibody description

MICALL2 Rabbit Polyclonal antibody. Positive WB detected in PC-3 cells. Observed molecular weight by Western-blot: 109 kDa

Preparation

This antibody was obtained by immunization of MICALL2 recombinant protein (Accession Number: NM_182924). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



PC-3 cells were subjected to SDS PAGE followed by western blot with Catalog No:112611(MICALL2 Antibody) at dilution of 1:1000