

USP37 antibody

Catalog Number: 116611

Product name

USP37 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

USP37 Rabbit Polyclonal antibody. Positive WB detected in HEK-293 cells, HeLa cells, HepG2 cells, human brain tissue, mouse brain tissue, mouse spleen tissue. Positive IP detected in mouse brain tissue. Observed molecular weight by Western-blot: 110 kDa

Preparation

This antibody was obtained by immunization of USP37 recombinant protein (Accession Number: BC112901). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:200-1:2000

Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:116611(USP37 antibody) at dilution of 1:500



IP Result of anti-USP37 (IP:Catalog No:116611, 4ug; Detection:Catalog No:116611 1:500) with mouse brain tissue lysate 4000ug.