

LPPR3 antibody



Catalog Number: 112311

Product name

LPPR3 antibody

Specificity

Human; other species not tested.

Antibody description

LPPR3 Rabbit Polyclonal antibody. Positive WB detected in K-562 cells, HL-60 cells. Observed molecular weight by Western-blot: 76 kDa

Preparation

This antibody was obtained by immunization of LPPR3 recombinant protein (Accession Number: BC017295). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



K-562 cells were subjected to SDS PAGE followed by western blot with Catalog No:112311(PRG2 Antibody) at dilution of 1:600