

Anti-Ooep antibody

Catalog Number: 176711

Product name

Anti-Ooep antibody

Specificity

Mouse

Antibody description

Rabbit polyclonal antibody to Ooep

Preparation

This antigen of this antibody was synthetic peptide within human ooop aa 52-96.

Formulation

Liquid, 1*PBS (pH7.4), 0.2% BSA, 40% Glycerol.
Preservative: 0.05% Sodium Azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ICC

Dilutions

ICC: 1:100

Validations




Fig1: Immunocytochemical staining of mouse two-cell embryos using anti-Ooep rabbit polyclonal antibody (red). (Green: F-actin, Blue: DAPI)