

VANGL1 antibody



Catalog Number: 116711

Product name

VANGL1 antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

VANGL1 Rabbit Polyclonal antibody. Positive WB detected in BxPC-3 cells. Observed molecular weight by Western-blot: 60-65 kDa

Preparation

This antibody was obtained by immunization of VANGL1 recombinant protein (Accession Number: NM_138959). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:1000

Validations



BxPC-3 cells were subjected to SDS PAGE followed by western blot with Catalog No:116711(VANGL1 antibody) at dilution of 1:300