

# CYP51A1 antibody



Catalog Number: 109711

## Product name

CYP51A1 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

CYP51A1 Rabbit Polyclonal antibody. Positive IP detected in mouse testis tissue. Positive WB detected in PC-3 cells, HL-60 cells, mouse kidney tissue, mouse lung tissue, mouse ovary tissue, mouse testis tissue. Observed molecular weight by Western-blot: 57kd

## Preparation

This antibody was obtained by immunization of CYP51A1 recombinant protein (Accession Number: NM\_000786). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:500-1:5000

## Validations



PC-3 cells were subjected to SDS PAGE followed by western blot with Catalog No:109711(CYP51A1 antibody) at dilution of 1:1000



IP Result of anti-CYP51A1 (IP:Catalog No:109711, 4ug; Detection:Catalog No:109711 1:1000) with mouse testis tissue lysate 2800ug.