

NANOS3 antibody



Catalog Number: 113021

Product name

NANOS3 antibody

Clonality

Polyclonal

Immunogen

[Human NANOS3 Recombinant protein \(GST tag\)](#)

Ig Type

Rabbit IgG

Specificity

Human, Mouse, Rat; other species not tested.

Applications

ELISA, WB

Antibody description

NANOS3 Rabbit Polyclonal antibody. Positive WB detected in mouse testis tissue, human kidney tissue, Transfected HEK-293 cells. Observed molecular weight by Western-blot: 26-30kd

Dilutions

Recommended Dilution:
WB: 1:200-1:1000

Preparation

This antibody was obtained by immunization of NANOS3 recombinant protein (Accession Number: NM_001098622). Purification method: Antigen affinity purified.


Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Validations



mouse testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:113021(NANOS3 antibody) at dilution of 1:300