

HS2ST1 Polyclonal Antibody



Catalog Number: 166021

Product name

Rabbit IgG

HS2ST1 Polyclonal Antibody

Applications

Specificity

WB

Human

Dilutions

Antibody description

WB 1:500 - 1:1000


Polyclonal antibody to HS2ST1

Validations

Preparation

HL-60

Antigen: A synthetic Peptide of human HS2ST1



Formulation

PBS with 0.02% sodium azide, pH7.3.

Storage

Store at 4°C. Avoid freeze / thaw cycles.

Western blot - HS2ST1 Polyclonal Antibody

Clonality

Polyclonal

Western blot analysis of extracts of HL-60 cells, using HS2ST1 antibody .Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution.Lysates/proteins: 25ug per lane.Blocking buffer: 3% nonfat dry milk in TBST.

Ig Type