

Anti-GNAT1 antibody

Catalog Number: 176021

Product name

Anti-GNAT1 antibody

Specificity

Human, Mouse, Rat

Antibody description

Rabbit polyclonal antibody to GNAT1

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human gnat1. aa range:71-120

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB

Dilutions

WB: 1:500-1:2000

Validations




Fig1.; Western Blot analysis of various cells using Gα t1 Polyclonal Antibody diluted at 1: 2000