

# Anti-Cdk4 Rabbit antibody

Catalog Number: 167021

## Product name

Anti-Cdk4 Rabbit antibody

WB: 1/1000

IHC: 1/50

## Specificity

Human, Mouse, Rat

ICC/IF: 1/50

IP: 1/20

## Antibody description

Cdk4 Rabbit polyclonal antibody

## Preparation

Antigen: Recombinant protein of human Cdk4

## Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

## Storage

Store at -20°C. Stable for 12 months from date of receipt.

## Clonality

Polyclonal

## Ig Type


Rabbit IgG

## Applications


WB, IHC-P, ICC/IF, IP

## Dilutions

## Validations



Western blot detection of Cdk4 in K562,C6,3T3,Hela cell lysates using Cdk4 Rabbit pAb(1:1000 diluted).Predicted band size:34kDa.Observed band size:34kDa.



Immunohistochemistry of Cdk4 in paraffin-embedded Human lung cancer tissue using Cdk4 Rabbit pAb at dilution 1/50