

# Anti-Hydroxyl-Histone H2A (Tyr39) Rabbit antibody



Catalog Number: 168021

## Product name

Monoclonal

Anti-Hydroxyl-Histone H2A (Tyr39) Rabbit antibody

## Ig Type

Rabbit IgG

## Specificity

Human, Mouse, Rat

## Applications

WB, IHC-P

## Antibody description

Hydroxyl-Histone H2A (Tyr39) Rabbit monoclonal antibody

## Dilutions

WB: 1/500-1/2000

IHC: 1/500-1/2000

## Preparation

Antigen: A synthesized peptide derived from human Histone H2A (hydroxyl Y39)

## Validations


## Formulation

Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality



Western blot detection of Hydroxyl-Histone H2A (Tyr39) in HeLa, A549, HL-60 cell lysates using Hydroxyl-Histone H2A (Tyr39) Rabbit mAb(1:1000 diluted). Observed band size: 14kDa.