

Anti-NEK7 Rabbit antibody



Catalog Number: 169021

Product name

Anti-NEK7 Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

NEK7 Rabbit monoclonal antibody

Preparation

Antigen: A synthetic peptide of human NEK7

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications


WB, IP

Dilutions

WB: 1/2000

IP: 1/20

Validations



Western blot detection of NEK7 in Rat Brain, C6 cell lysates using NEK7 Rabbit mAb(1:1000 diluted). Predicted band size: 35kDa. Observed band size: 35kDa.