

# Human UBE2A/RAD6A (GST & His tag) recombinant protein



Catalog Number: 519021

## General Information

### Gene Name Synonym

Ubiquitin-Conjugating Enzyme E2 A, RAD6 Homolog A, HR6A, hHR6A, Ubiquitin Carrier Protein A, Ubiquitin-Protein Ligase A, UBE2A, RAD6A

### Protein Construction

A DNA sequence encoding the Human UBE2A/RAD6A Met1-Cys152 is expressed with a GST, 6His tag at the N-terminus.

### Organism

Human

### Expression Host

E. coli

## QC Testing

### Purity

>95% as determined by SDS-PAGE.

### Endotoxin

<1.0 EU per µg protein as determined by the LAL

method.

### Stability

Store for up to 6 months at -20°C

### Formulation


Supplied as a 0.2 µm filtered solution of 50mM HEPES, 150mM NaCl, 2mM DTT, 10% Glycerol, pH 7.5.

## Usage Guide

### Storage

Store at <-20°C, minimize freeze-thaw cycles.

### SDS-PAGE



Human UBE2A/RAD6A (GST & His tag)  
recombinant protein