

HLA-DPB1 antibody



Catalog Number: 111421

Product name

HLA-DPB1 antibody

Specificity

Human; other species not tested.

Antibody description

HLA-DPB1 Rabbit Polyclonal antibody. Positive IHC detected in human spleen tissue, human tonsillitis tissue.

Preparation

This antibody was obtained by immunization of HLA-DPB1 recombinant protein (Accession Number: NM_002121). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IHC

Dilutions

Recommended Dilution:

IHC: 1:20-1:200

Validations



Immunohistochemistry of paraffin-embedded human spleen tissue slide using Catalog No:111421(HLA-DPB1 Antibody) at dilution of 1:50 (under 40x lens)