

GSTM3 Polyclonal Antibody



Catalog Number: 162421

Product name

GSTM3 Polyclonal Antibody

Specificity

Human

Antibody description

Polyclonal antibody to GSTM3

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-225 of human GSTM3 (NP_000840.2).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

IF

Dilutions

IF 1:50 - 1:200

Validations



Immunofluorescence - GSTM3 Polyclonal Antibody

Immunofluorescence analysis of U2OS cells using GSTM3 antibody at dilution of 1:100. Blue: DAPI for nuclear staining.