

# MCRS1 antibody



Catalog Number: 107421

## Product name

MCRS1 antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

MCRS1 Mouse Monoclonal antibody. Positive WB detected in Skov-3 cells. Observed molecular weight by Western-blot: 55+110KD

## Preparation

This antibody was obtained by immunization of MCRS1 recombinant protein (Accession Number: NM\_006337). Purification method: Protein A purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Monoclonal

## Ig Type

Mouse IgG2a

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



SKOV-3 cells were subjected to SDS PAGE followed by western blot with Catalog No:107421(MCRS1 antibody) at dilution of 1:1000