

Anti-Cytokeratin 19 (1H6) Mouse antibody



Catalog Number: 168421

Product name

Mouse IgG1

Anti-Cytokeratin 19 (1H6) Mouse antibody

Applications

Specificity

WB, IHC-P

Human

Dilutions

Antibody description

WB: 1/500 - 1/2000

Cytokeratin 19 (1H6) Mouse monoclonal antibody

IHC: 1/200 - 1/1000

Preparation

ELISA: 1/10000

Antigen: Purified recombinant fragment of human KRT19 expressed in E. Coli.

Validations

Formulation

Ascitic fluid containing 0.03% sodium azide.


Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type



Western blot detection of Cytokeratin 19 in CaCo2, HepG2 and T-47D cell lysates using Cytokeratin 19 mouse mAb (dilution 1:500). Predicted band size: 41 kDa. Observed band size: 41 kDa.