

GADD45G antibody



Catalog Number: 110821

Product name

GADD45G antibody

Specificity

Human; other species not tested.

Antibody description

GADD45G Rabbit Polyclonal antibody. Positive WB detected in UV treated HeLa cells, human placenta tissue. Observed molecular weight by Western-blot: 17kd

Preparation

This antibody was obtained by immunization of GADD45G recombinant protein (Accession Number: NM_006705). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:1000

Validations



UV treated HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:110821(GADD45G Antibody) at dilution of 1:300